



## Pharmenta® AptiPort® Sampling Valve 1.5TC-W

### Notes:

Mounts to flush with tank wall through a 1.5" sanitary tri-clamp tank port.

- Fits into and passively drains from any 1.5" sanitary tri-clamp tank port.
- Mounts flush with inside tank wall on 1.5" tri-clamp ports up to 1.4" long.

**Pharmenta**  
Princeton, New Jersey

THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF NL Technologies Ltd. ANY REPRODUCTION IN PART OR WHOLE WITHOUT THE WRITTEN PERMISSION OF NL Technologies Ltd. IS PROHIBITED.

**APTIPORT VALVE,  
1 INCH, TC - REAR**

SIZE <b>A</b>	DWG. NO. <b>100196</b>	REV. <b>A</b>
SHT 1 OF 2		